

# SCALEX™ 850

## SELF-CONTAINED ULTRASONIC SCALER



SCALEX™ 850

- Self-contained solution tank allows for convenient portability.
- One touch selection between medicament tank and external water source.
- 25k/30k scaling strokes per second, solid state power unit.
- Innovative water-cooled circuitry designed for greater stability during operation.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit is available in 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001  
ISO 13485



**DENTAMERICA®**

# SCALEX™ 850

## SELF-CONTAINED ULTRASONIC SCALER

### ※ULTRASONIC INSERTS

**Reorder No.80005**


Universal Insert 25K  
(External Water)

**Reorder No.83005**


Universal Insert 30K  
(External Water)

**Reorder No.80006**


Universal Insert 25K  
(Internal Water)

**Reorder No.83006**


Universal Insert 30K  
(Internal Water)

**Reorder No.80007**


Chisel Type Insert 25K  
(External Water)

**Reorder No.83007**


Chisel Type Insert 30K  
(External Water)

**Reorder No.80008**

Chisel Type Insert 25K  
(Internal Water)

**Reorder No.83008**

Chisel Type Insert 30K  
(Internal Water)

**Reorder No.80011**

Endodontic Insert 25K  
(Canalsonic Tip)

**Reorder No.88060****Foot control**

Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

**SPECIFICATIONS:**

- \* FREQUENCY: 25,000 or 30,000 Hz/sec.
- \* DIMENSIONS: D260xW140xH210 mm
- \* VOLUME: 1,000 ml
- \* WEIGHT: 3.3 kg
- \* INPUT: 115V/60Hz or 230V/50Hz
- \* OUTPUT: AC 24V 3A

**DENTAMERICA®**

18688 E.San Jose Ave.  
City of Industry, CA 91748 U.S.A.  
Phone: (626) 912-1388  
Fax: (626) 913-0510  
www.dentamerica.com  
www.dentamerica.asia

**Reorder No.85003****Self-contained solution tank**

Easily change the solutions according to different applications including periodontics, endodontics, and prophylaxis.

DISTRIBUTOR