

SCALEX™ 830

ULTRASONIC SCALER



SCALEX™
830

- 30,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation.
- No power loss at maximum pressure.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit for 25K scaling strokes is available.
- Unit is available in 110/120V 50/60Hz, or 220/240V 50Hz operation.

ISO 9001
ISO 13485



GMP

DENTAMERICA®

SCALEX™ 830

ULTRASONIC SCALER

※ULTRASONIC INSERTS

Reorder No.80005

Universal Insert 25K
(External Water)

**Reorder No.83005**

Universal Insert 30K
(External Water)

**Reorder No.80006**

Universal Insert 25K
(Internal Water)

**Reorder No.83006**

Universal Insert 30K
(Internal Water)

**Reorder No.80007**

Chisel Type Insert 25K
(External Water)

**Reorder No.83007**

Chisel Type Insert 30K
(External Water)

**Reorder No.80008**

Chisel Type Insert 25K
(Internal Water)

**Reorder No.83008**

Chisel Type Insert 30K
(Internal Water)

**Reorder No.80011**

Endodontic Insert 25K
(Canalsonic Tip)

**SPECIFICATIONS :**

- INPUT : 115V 60Hz or 230V 50Hz 3A
- OUTPUT : AC 24V 2.4A
- FREQUENCY : 25,000 or 30,000 Hz/sec.
- DIMENSION : W140xD140xH50mm
- WEIGHT : 1000g

**Reorder
No.88060****Foot control**

Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

DENTAMERICA®

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com

DISTRIBUTOR